

# CLINIVIEW™ OrthoTrace



Perfect orthodontic tool for cephalometric analysis and treatment planning



## CLINIVIEW™ OrthoTrace

Flexible and user-friendly software for orthodontists and maxillofacial surgeons

- Orthodontic tracing
- Cephalometric analysis (lateral and PA)
  - More than 30 standard analyses available
  - Customized analysis
- Superimposition
- Treatment planning
- Soft tissue prediction
  - Including profile photo prediction
- Cast model analysis



# CLINIVIEW™ OrthoTrace

## Perfect orthodontic tool for cephalometric analysis and treatment planning

- Orthodontic tracing • Cephalometric analysis • Superimposition • Treatment planning • Soft tissue prediction

### TECHNICAL SPECIFICATIONS

#### Basic features

Supports Windows 8/7/Vista operating systems

Direct integration with CLINIVIEW™ 10.2.3 (or later) and connection also to patient management systems or PACS

Import of digital images (JPEG, TIFF, BMP and DICOM)

Image calibration

Image functions (zoom, brightness/contrast, rotate, mirror, flip, invert)

#### Cephalometry

More than 30 standard analyses available

Easy-to-use analysis editor for customized analysis

Interactive measurements

Cast model analysis

#### Tracing

Place markers (landmarks)

Marker guide

Automatic placement of incisors

Automatic positioning of the cursor to the landmark

Place/Draw teeth and hard tissue

Place/Draw soft tissue profile line

#### Superimposition

Superimpose tracings

Superimpose radiographs (can be done based on stable anatomical structures)

Superimpose lateral radiograph with profile photo

#### Output

Printing

Reports

Copy/Export images and data

#### VTO – Visual Treatment Objective

Planned movement of hard tissue and teeth

Soft tissue profile prediction

Profile photo prediction

#### Miscellaneous

Multiple language support

Documentation in pdf



© 2014 Instrumentarium Dental, 211149-1 English  
CLINIVIEW™ is a common law trademark of Instrumentarium Dental, PaloDEX Group Oy.

#### Instrumentarium Dental

Nahkelantie 160, P.O. Box 20 FI-04301 Tuusula, Finland  
Tel. +358 10 270 2000  
www.instrumentariumdental.com, info@instrudental.com



INSTRUMENTARIUM